## ENDOBUTTON® Direct Fixation Device

### ENDOBUTTON Direct-Fixationsinstrument

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#### English

#### English

## ENDOBUTTON® Direct Fixation Device

### Device Description

The Smith & Nephew ENDOBUTTON® Direct Fixation Device provides the orthopedic surgeon a means of accurate cortical fixation in reconstructive surgery. The fixation device allows for endoscopic ligament reconstruction without an ancillary lateral incision. It provides an alternative fixation method for soft tissue grafts, such as semitendinosus, which do not include bone block. The fixation device, available for graft sizes 5–10 mm, is made of titanium alloy, which meets ASTM standard F336. The fixation device is available with and without sutures. The preloaded suture consists of a #5 polyethylene braided suture (lead suture) and a #5 polyester braided suture (flipping suture).

The fixation device is provided sterile, for single use only.

### Indications for Use

The Smith & Nephew ENDOBUTTON Direct Fixation Device is used for fixation of tendons and ligaments during orthopedic reconstruction procedures such as anterior cruciate ligament (ACL) and posterior cruciate ligament (PCL).

### Contraindications

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Insufficient quantity or quality of bone.
- Blood supply or previous infections which may tend to retard healing.
- Active infection.
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.

### Warnings

-  Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.
- Contents are sterile unless package is opened or damaged. DO NOT RESTERILIZE. For single use only. Discard any open, unused product. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- The ENDOBUTTON Direct Fixation Device has not been evaluated for safety and migration in the magnetic resonance (MR) environment. The ENDOBUTTON Direct Fixation Device has not been tested for heating or migration in the MR environment.

### Precautions

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- Hazards associated with reuse of this device include, but are not limited to, patient infection and/or device malfunction.
- Prior to use, inspect the product package for signs of damage or tampering. If damaged, do not use.
- Prior to use, examine the device for possible damage to assure proper functioning. If damaged, do not use.
- The use of metallic surgical implants provides the orthopedic surgeon with a means of accurate fixation and helps generally in the management of fractures and reconstructive surgery. These implants are intended as aids to normal healing, but are not intended to replace normal body structures or bear the weight of the body in the presence of incomplete bone healing<sup>†</sup>.
- Postoperative care is important. A patient should be instructed on the limitations of the implant and should be cautioned regarding weight bearing and body stresses on the appliance prior to secure bone healing.
- Additional precautions include those applicable to arthroscopy. In general, careful attention must be paid to asepsis and avoidance of anatomical hazards.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

### Adverse Reactions

Complications are those seen with any method of internal fixation.

### Instructions for Use

For details refer to the Smith & Nephew ACL Reconstruction with the ACUFEX® Director Drill Guide and ENDOBUTTON CL Fixation System Surgical Technique (REF 10600026).

- Verify graft size (between 5 and 10 mm).
- Prepare the tibial tunnel using standard technique.
- Drill the femoral tunnel (5–10 mm diameter) through the cortex and measure its length.
- Choose an appropriately sized fixation device and load the graft as indicated for each device.
- ENDOBUTTON Direct Fixation Device (without sutures): Attach a #5 braided polyester suture (lead suture) to the ENDOBUTTON Direct Fixation Device in the hole on the closed end.
- Attach a #5 flipping suture to the ENDOBUTTON Direct Fixation Device through the two holes at the open end (Figure 1).


**Figure 1**  
Alternatively, secure the trailing suture to the graft to keep it centered to the button (Figures 2 and 3).


**Figure 2**



**Figure 3**  
Using a marker, make a line on the graft at the insertion point of the femoral tunnel. Then place a second line (according to chart below) distal to the first line to indicate the flipping distance (Figure 4).

| Tunnel Size | Marking |
|-------------|---------|
| 5 mm        | 8 mm    |
| 6 mm        | 9 mm    |
| 7 mm        | 10 mm   |
| 8 mm        | 11 mm   |
| 9 mm        | 11 mm   |
| 10 mm       | 12 mm   |

- Pull the graft into place using the standard ENDOBUTTON device technique.  
**Note:** To prevent pulling the device into the soft tissue, apply sufficient counter tension to the graft on the tibial side.  
**Note:** Reference the second line on the graft to ensure the graft is not pulled out soft tissue.
- Rotate the fixation device immediately external to the femur, locking in place as the graft retreats. Verify that the first mark on graft is properly positioned at the femoral tunnel entrance inside the joint.  
X-ray or fluoroscopy confirms the position of the fixation device on the anterolateral femoral cortex.

### Warranty

For single use only. This product is warranted to be free from defects in material and workmanship. Do not reuse.

### For Further Information

If further information on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

*\* Refer to product-specific packaging to determine actual sterilization method.*

<sup>†</sup>Trademark of Smith & Nephew, registered U.S. Patent & Trademark Office.  
Covered by U.S. patent numbers 5,306,301, 5,645,588, 5,769,894, and 7,530,990. Other patents pending.  
©2011 Smith & Nephew, Inc. All rights reserved.

#### English

#### English

#### English

## ENDOBUTTON® Direct-Fixationsinstrument

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\*Varumärke som tillhör Smith & Nephew, och har registrerats vid U.S. Patent and Trademark Office (USA:s Patent- och registreringsbyrå). Omtalats av de amerikanska patentnumren 5.306.301, 5.645.588, 5.769.894, och 7.530.990. Andra patentansökningar har inlämnats ©2011 Smith & Nephew, Inc. Med ensamrätt.